

ELECTRIC GYNECOLOGICAL TABLE HBD34-EG



ELECTRIC GYNECOLOGICAL TABLE HBD34-EG

The Delivery table is a flexible designed for suitable in conducting parturient delivery, gynecological surgery, operative abortion, diagnosis and examinations etc. by medical unit. It has an array of user-friendly features. The table has been designed keeping in mind the safety and comfort of the patient to support throughout the stages of delivery. Used in Gyneacological surgery, Operative abortion, Gyneacological examination, Labor, delivery, postpartum recovery period, Parturient delivery, Diagnosis, Ultrasonic examination, Gynecological abdominal operation, Vulva operation. Also known as Gyneacological table, Gyneacological bed, Examination table, Examination bed, Delivery bed, Gynaecological Delivery Bed.

HBD34-EG ELECTRIC GYNECOLOGICAL TABLE

Various actions are controlled by Imported brand motor

Hand-held controller and foot controller, mutually independent and is easily for operation

Leg support frame can fold up and down electric

Unique profile design which is beautiful and elegant

Whole stainless steel hidden auxiliary table with flexible extension and retraction

Equipped with hidden sewage basin

The table adopts high density memory sponge and is once foam molding, with no seam

Features with one-key reset function

Super mute foot wheel and equipped with brake device



SPECIFICATIONS

Model	HBD34-EG
Back-rest lifting angle	$\geq 40^\circ$
Foot-rest lifting angle	$\geq 30^\circ$
Trendelenburg	$\geq 5^\circ$
Reverse trendelenburg	$\geq 12^\circ$
Working area	$(350 \times 480) \pm 20$
Dimension	1270Lx550Wx700~1000 mm
Voltage	220V \pm 22V, 50Hz \pm 1Hz



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com