

# ELECTRIC OBSTERIC BED HBD35-EO



# ELECTRIC OBSTERIC BED HBD35-EO

The Delivery table is a flexible designed for suitable in conducting parturient delivery, gynecological surgery, operative abortion, diagnosis and examinations etc. by medical unit. It has an array of user-friendly features. The table has been designed keeping in mind the safety and comfort of the patient to support throughout the stages of delivery.

Used in Gyneacological surgery, Operative abortion, Gyneacological examination, Labor, delivery, postpartum recovery period, Parturient delivery, Diagnosis, Ultrasonic examination, Gynecological abdominal operation, Vulva operation. Also known as Gyneacological table, Gyneacological bed, Examiniation table, Examiniation bed, Delivery bed, Gynaecological Delivery Bed.

## HBD35-EO ELECTRIC OBSTERIC BED

Imported LINAK motor various actions are controlled by remote controller

Super mute big foot wheel which is easily moved and equipped with center control breaking system.

Equipped with stainless steel sewage basin

The mattress adopts high density sponge and manufactured by high quality artificial leather without seam

With various complete accessories

Standard Accessories: A pair of leg support frame, a stainless steel sewage basin, a stand for transfusion, a remote control, a power wire



## SPECIFICATIONS

Model	HBD35-EO
Back-rest lifting angle	$\geq 80^\circ$
Foot-rest lifting angle	$\geq 45^\circ$
Trendelenburg	$\geq 5^\circ$
Reverse trendelenburg	$\geq 15^\circ$
Dimension	2010Lx780Wx565-935 mm



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)