

MANUAL MULTI- FUNCTION GYNEACOLOGICAL TABLE HBD36-MM



MANUAL MULTI-FUNCTION GYNEACOLOGICAL TABLE HBD36-MM

The Delivery table is a flexible designed for suitable in conducting parturient delivery, gynecological surgery, operative abortion, diagnosis and examinations etc. by medical unit. It has an array of user-friendly features. The table has been designed keeping in mind the safety and comfort of the patient to support throughout the stages of delivery.

Used in Gyneacological surgery, Operative abortion, Gyneacological examination, Labor, delivery, postpartum recovery period, Parturient delivery, Diagnosis, Ultrasonic examination, Gynecological abdominal operation, Vulva operation.

Also known as Gyneacological table, Gyneacological bed, Examination table, Examination bed, Delivery bed, Gynaecological Delivery Bed.

HBD36-MM MANUAL MULTI-FUNCTION GYNEACOLOGICAL TABLE

Material: Steel Powder Coated S.S, Color: Blue



SPECIFICATIONS

Model	HBD36-MM
Function	Multi
Back-rest lifting angle	$\geq 20^\circ$
Back section up adjusted	$\geq 60^\circ$
Seat board upwards	$\geq 8^\circ$
Seat board downwards	$\geq 10^\circ$
Load Bearing	250 kg
Height adjustment	≤ 650 mm
Tabletop Dimension	1240Lx600W mm
Upright Dimension	500Lx600W mm
Dimension	1750Lx600Wx650~900 mm



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: info@biozef.com | Website: biozef.com