

STAINLESS STEEL CHILDREN BED HBD41-SS



STAINLESS STEEL CHILDREN BED HBD41-SS

Children medical Bed built to the most current standards while also meeting the clinician's needs which allows easy care and disinfection. The bed is equipped with Trendelenburg and reverse Trendelenburg positions. Simple and modern design for safety and protection of the child.

Used in Medical unit, Nursing.

Also known as Children Hospital Bed, Medical bed, Hospital bed, Medical child bed, Pediatric bed, Medical pediatric bed, Hospital pediatric bed, Children medical bed, Children pediatric hospital bed.

HBD41-SS STAINLESS STEEL CHILDREN BED

Stainless steel framework

Stainless steel side rail, drop type

Noiseless castors

With mattress

With debris basket



SPECIFICATIONS

Model	HBD41-SS
IV pole	Yes
Dimension	750x450x610-1030 mm



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com