

# CHILDREN BED WITH TWO FUNCTION HBD47-02



# CHILDREN BED WITH TWO FUNCTION HBD47-02

Children medical Bed built to the most current standards while also meeting the clinician's needs which allows easy care and disinfection. The bed is equipped with Trendelenburg and reverse Trendelenburg positions. Simple and modern design for safety and protection of the child.

Used in Medical unit, Nursing.

Also known as Children Hospital Bed, Medical bed, Hospital bed, Medical child bed, Pediatric bed, Medical pediatric bed, Hospital pediatric bed, Children medical bed, Children pediatric hospital bed.

## HBD47-02 CHILDREN BED WITH TWO FUNCTION



### SPECIFICATIONS

|                             |   |
|-----------------------------|---|
| Model                       | HBD47-02                                    |
| Load Bearing                | 150 kg                                      |
| Dimension                   | 1770x900 mm                                 |
| Packing Dimension           | wooden carton: 1880x980x1450 mm (assembled) |
| carton                      |   |
| 1880x1020x310 mm (detached) | Net Weight                                  |



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)