

STAINLESS STEEL PATIENT TROLLEY HBD81-SS



STAINLESS STEEL PATIENT TROLLEY HBD81-SS

Labstac Stretcher is an economical choice designed with ease of use and versatility in mind. Adjustable heights for working in extreme operative conditions. A stretcher that allows the operator to assist the patients wherever they may be in a fast and secure way. Exceptional geometrical design provides you with confidence in the durability of our stretcher fleet.

Used in Patient transfer, Patient trolley.

Also known as Medical stretcher, Ambulance stretcher, Emergency Stretcher, Patient transfer stretcher.

HBD81-SS STAINLESS STEEL PATIENT TROLLEY

Stainless steel frame

5" casters with brake

Stainless steel guardrail can be putted down

Separated into two parts: stretcher and trolley



SPECIFICATIONS

Model	HBD81-SS
Dimension	1960x600x750 mm



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com