

# HYDRAULIC TRANSPORTATION TROLLEY HBD82-HT



# HYDRAULIC TRANSPORTATION TROLLEY HBD82-HT

Labstac Stretcher is an economical choice designed with ease of use and versatility in mind. Adjustable heights for working in extreme operative conditions. A stretcher that allows the operator to assist the patients wherever they may be in a fast and secure way. Exceptional geometrical design provides you with confidence in the durability of our stretcher fleet.

Used in Patient transfer, Patient trolley.

Also known as Medical stretcher, Ambulance stretcher, Emergency Stretcher, Patient transfer stretcher.

## HBD82-HT HYDRAULIC TRANSPORTATION TROLLEY

Material: Galvanised steel tube and surface coated by electrostatic spray

The back board is X-ray translucent, with X-ray cassette

Hydraulic pump imported from USA Actuant

Foot pedal control at both sides of the trolley

Advanced centrally controlled brake system with 200 mm dustproof castors

The bed can be raised and lowered, tilted forward and backward

The bottom cover is made of PP material, with oxygen cylinder holder

With stainless steel foldable railing

With 5th directional castors

with infusion stand, drainage hook

With foam mattress PU material 5 cm thickness



## SPECIFICATIONS

Model	HBD82-HT
Dimension	2000x650x630-930 mm



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)