

# STAINLESS STEEL PATIENT TROLLEY HBD83-SS



# STAINLESS STEEL PATIENT TROLLEY HBD83-SS

Labstac Stretcher is an economical choice designed with ease of use and versatility in mind. Adjustable heights for working in extreme operative conditions. A stretcher that allows the operator to assist the patients wherever they may be in a fast and secure way. Exceptional geometrical design provides you with confidence in the durability of our stretcher fleet.

Used in Patient transfer, Patient trolley.

Also known as Medical stretcher, Ambulance stretcher, Emergency Stretcher, Patient transfer stretcher.

## HBD83-SS STAINLESS STEEL PATIENT TROLLEY

Manual head and foot restricted

Stainless steel frame

Aluminium alloy railing

With mattress



## SPECIFICATIONS

Model	HBD83-SS
Dimension	1930x750x620-880 mm



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)