

OPERATION CONNECTING TROLLEY HBD84-OC



OPERATION CONNECTING TROLLEY HBD84-OC

Labstac Stretcher is an economical choice designed with ease of use and versatility in mind. Adjustable heights for working in extreme operative conditions. A stretcher that allows the operator to assist the patients wherever they may be in a fast and secure way. Exceptional geometrical design provides you with confidence in the durability of our stretcher fleet.

Used in Patient transfer, Patient trolley.

Also known as Medical stretcher, Ambulance stretcher, Emergency Stretcher, Patient transfer stretcher.

HBD84-OC OPERATION CONNECTING TROLLEY

Material: made of quality Aluminium alloy. Cart surface and side rails are shaped by once injecting molding of PP engineering plastics

The guideways for linking two small separate carts are made of quality Aluminium alloy of high strength and thus enable the bed to move steadily and reliably

Two separate small carts can be raised and lowered for adjustment. Only when the two small carts are adjusted to the same level, the bed rest can be slid smoothly

As the bed slides to any separate cart, it can be automatically locked and equipped with safety device

An advanced centrally controlled brake system is used for both separate small carts in order to make them smoother and more stable when linking up and running



SPECIFICATIONS

| | |
|-----------|---------------------|
| Model | HBD84-OC |
| Dimension | 3650x640x640-970 mm |



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com