

# HOSPITAL TRANSFER STRETCHER TROLLEY HBD8C-00



# HOSPITAL TRANSFER STRETCHER TROLLEY

## HBD8C-00

Labstac Stretcher is an economical choice designed with ease of use and versatility in mind. Adjustable heights for working in extreme operative conditions. A stretcher that allows the operator to assist the patients wherever they may be in a fast and secure way. Exceptional geometrical design provides you with confidence in the durability of our stretcher fleet.

Used in Patient transfer, Patient trolley.

Also known as Medical stretcher, Ambulance stretcher, Emergency Stretcher, Patient transfer stretcher.

## HBD8C-00 HOSPITAL TRANSFER STRETCHER TROLLEY

It is a type of manual transfer stretcher. Turn by deasil, the stretcher will be up, and turn by widdershins, the stretcher will be down;

Aluminium-Alloy protective rails, can be folded and hidden under the surface on no use;

Central-lock castors, diameter 150mm. Matched with direction castor, help the operator controlling the direction easily;

Rotate the crand and handle on bedboard, patient on the mattress can be moved parallelly to another bed or flat table



## SPECIFICATIONS

Model	HBD8C-00
Height adjustment	920-610 mm
Net Weight	89 kg



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)