

MANUAL TRANSFER STRETCHER HBD8G- MS1



MANUAL TRANSFER STRETCHER HBD8G-MS1

Labstac Stretcher is an economical choice designed with ease of use and versatility in mind. Adjustable heights for working in extreme operative conditions. A stretcher that allows the operator to assist the patients wherever they may be in a fast and secure way. Exceptional geometrical design provides you with confidence in the durability of our stretcher fleet.

Used in Patient transfer, Patient trolley.

Also known as Medical stretcher, Ambulance stretcher, Emergency Stretcher, Patient transfer stretcher.

HBD8G-MS1 MANUAL TRANSFER STRETCHER

Single crank for height adjustment Turn by clockwise, stretcher will be up; Turn by anti- clockwise, stretcher will be down

Bright arrow symbols direct the user how to operate clearly

Central-lock castors with steering, diameter 150mm

Concave plastic chassis cover, for storage use

Rotary rollover PP protective railing is blow molding by one-shot plastic forming, 3 different positions of the side rail supply different use

PP bed board produced by blowing molding, easy for clean

Standard accessories: Mattress, I.V. pole



SPECIFICATIONS

Model	HBD8G-MS1
Back-rest lifting angle	0-70°±5°
Height adjustment	795- 510 mm (from floor)
IV pole	Yes
Dimension	1960Lx680W mm
Net Weight	67 kg



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com