

CONNECTING TRANSFER STRETCHER HBD8H-CS



CONNECTING TRANSFER STRETCHER HBD8H-CS

Labstac Stretcher is an economical choice designed with ease of use and versatility in mind. Adjustable heights for working in extreme operative conditions. A stretcher that allows the operator to assist the patients wherever they may be in a fast and secure way. Exceptional geometrical design provides you with confidence in the durability of our stretcher fleet.

Used in Patient transfer, Patient trolley.

Also known as Medical stretcher, Ambulance stretcher, Emergency Stretcher, Patient transfer stretcher.

HBD8H-CS CONNECTING TRANSFER STRETCHER

It is a manual Connecting Transfer Stretcher, can be used to defend cross-infection

Cart A inside the operation room, Cart B outside, adjusting them to the same height, then push Cart B to Cart A, two carts will be connected, following the orbit, bedboard C can be pushed from one cart to another, and VICE VERSA

Bright arrow symbols direct the users how to operate clearly

Standard accessories: Mattress, I.V. pole



SPECIFICATIONS

Model	HBD8H-CS
Back-rest lifting angle	0-60°
Height adjustment	915-660 mm
IV pole	Yes
Net Weight	139 kg



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com