

HYDRAULIC TRANSFER STRETCHER HBD8I-HS



HYDRAULIC TRANSFER STRETCHER HBD8I-HS

Labstac Stretcher is an economical choice designed with ease of use and versatility in mind. Adjustable heights for working in extreme operative conditions. A stretcher that allows the operator to assist the patients wherever they may be in a fast and secure way. Exceptional geometrical design provides you with confidence in the durability of our stretcher fleet.

Used in Patient transfer, Patient trolley.

Also known as Medical stretcher, Ambulance stretcher, Emergency Stretcher, Patient transfer stretcher.

HBD8I-HS HYDRAULIC TRANSFER STRETCHER

It is a type of hydraulic transfer stretcher

Five functions available by pressing the foot-pedals

Bight arrow symbols direct the users how to operate clearly

Central-lock castors, diameter 150mm. Matched with direction castor under the chassis, which help users easily control the stretcher direction

Standard accessories: Mattress, I.V.Pole



SPECIFICATIONS

Model	HBD8I-HS
Back-rest lifting angle	0-60°
Trendelenburg	0-13°
Reverse trendelenburg	0-13°
Height adjustment	815-580 mm
IV pole	Yes
Net Weight	82 kg



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com