

HYDRAULIC TRANSFER STRETCHER HBD8J- HS1



HYDRAULIC TRANSFER STRETCHER HBD8J-HS1

Labstac Stretcher is an economical choice designed with ease of use and versatility in mind. Adjustable heights for working in extreme operative conditions. A stretcher that allows the operator to assist the patients wherever they may be in a fast and secure way. Exceptional geometrical design provides you with confidence in the durability of our stretcher fleet.

Used in Patient transfer, Patient trolley.

Also known as Medical stretcher, Ambulance stretcher, Emergency Stretcher, Patient transfer stretcher.

HBD8J-HS1 HYDRAULIC TRANSFER STRETCHER

It is a type of hydraulic transfer stretcher, 5 functions available by hydraulic and gas spring and crank system together

At the head-end, there is a pair of push- handles, can be folded up and down

Bright arrow symbols direct the users how to operate clearly

Central-lock castors, diameter 150mm. Matched with direction castor, help the operator controlling the direction easily, Bed board for X-Ray Examination



SPECIFICATIONS

Model	HBD8J-HS1
Back-rest lifting angle	0-60°
Foot-rest lifting angle	0-40°
Trendelenburg	0-13°
Reverse trendelenburg	0-13°
Height adjustment	835-600 mm
Net Weight	105 kg



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com