

MANUAL WHEELCHAIR HBD91- MW0



MANUAL WHEELCHAIR HBD91-MW0

Wheelchair is mainly used as a manually operated or power-driven device designed with double cross-shaped brackets, which enhance the safety of the structure. In the meantime, adjustable armrests, quick disassembly of footrest, positioning of Controller and other functions reflect the user-friendly design, and provide great convenience for the user.

Used in Patient Indoor-outdoor locomotion, Patient transfer.

Also known as Patient transfer wheelchair, Transport wheelchairs.

HBD91-MW0 MANUAL WHEELCHAIR

It is a type of steady and simple Manual Wheelchair



SPECIFICATIONS

Model	HBD91-MW0
Wheelchair length	970 mm
Opened wheelchair width	660 mm
Wheelchair Height	920 mm
Closed wheelchair width	350 mm
Seat Width	430 mm
Front wheel	6" inch
Rear wheel	24" inch
Seat Height	470 mm
Seat depth	400 mm
Backrest Height	450 mm
Armrest Height	210 mm
Distance between armrest	425 mm
Height between footrests and ground	120 mm
Net Weight	12.7 Kg



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: info@biozef.com | Website: biozef.com