

# MANUAL WHEELCHAIR HBD91- MW2



# MANUAL WHEELCHAIR HBD91-MW2

Wheelchair is mainly used as a manually operated or power-driven device designed with double cross-shaped brackets, which enhance the safety of the structure. In the meantime, adjustable armrests, quick disassembly of footrest, positioning of Controller and other functions reflect the user-friendly design, and provide great convenience for the user.

Used in Patient Indoor-outdoor locomotion, Patient transfer.

Also known as Patient transfer wheelchair, Transport wheelchairs.

## HBD91-MW2 MANUAL WHEELCHAIR

It is a type of steady and simple Manual Wheelchair



### SPECIFICATIONS

Model	HBD91-MW2
Wheelchair length	1020 mm
Opened wheelchair width	630 mm
Wheelchair Height	890 mm
Closed wheelchair width	345 mm
Seat Width	430 mm
Front wheel	6" inch
Rear wheel	24" inch
Seat Height	460 mm
Seat depth	400 mm
Backrest Height	450 mm
Armrest Height	210 mm
Distance between armrest	445 mm
Height between footrests and ground	80 mm
Net Weight	13.28 Kg





**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)