

# MANUAL WHEELCHAIR HBD91- MW4



# MANUAL WHEELCHAIR HBD91-MW4

Wheelchair is mainly used as a manually operated or power-driven device designed with double cross-shaped brackets, which enhance the safety of the structure. In the meantime, adjustable armrests, quick disassembly of footrest, positioning of Controller and other functions reflect the user-friendly design, and provide great convenience for the user.

Used in Patient Indoor-outdoor locomotion, Patient transfer.

Also known as Patient transfer wheelchair, Transport wheelchairs.

## HBD91-MW4 MANUAL WHEELCHAIR

It is a type of steady and simple Manual Wheelchair



### SPECIFICATIONS

Model	HBD91-MW4
Wheelchair length	855 mm
Opened wheelchair width	535 mm
Wheelchair Height	855 mm
Closed wheelchair width	250 mm
Seat Width	450 mm
Front wheel	6" inch
Rear wheel	14" inch
Seat Height	480 mm
Seat depth	430 mm
Backrest Height	350 mm
Armrest Height	190 mm
Distance between armrest	440 mm
Height between footrests and ground	60 mm
Net Weight	13.75 Kg





**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)