

MANUAL WHEELCHAIR HBD91- MW7



MANUAL WHEELCHAIR HBD91-MW7

Wheelchair is mainly used as a manually operated or power-driven device designed with double cross-shaped brackets, which enhance the safety of the structure. In the meantime, adjustable armrests, quick disassembly of footrest, positioning of Controller and other functions reflect the user-friendly design, and provide great convenience for the user.

Used in Patient Indoor-outdoor locomotion, Patient transfer.

Also known as Patient transfer wheelchair, Transport wheelchairs.

HBD91-MW7 MANUAL WHEELCHAIR

It is a type of steady and simple Manual Wheelchair



SPECIFICATIONS

Model	HBD91-MW7
Wheelchair length	830 mm
Opened wheelchair width	580 mm
Wheelchair Height	780 mm
Closed wheelchair width	350 mm
Seat Width	360 mm
Front wheel	4" inch
Rear wheel	24" inch
Seat Height	470 mm
Seat depth	380 mm
Backrest Height	330 mm
Armrest Height	215 mm
Distance between armrest	415 mm
Height between footrests and ground	100 mm
Net Weight	13 Kg



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: info@biozef.com | Website: biozef.com