

# MANUAL WHEELCHAIR HBD91- MW8



# MANUAL WHEELCHAIR HBD91-MW8

Wheelchair is mainly used as a manually operated or power-driven device designed with double cross-shaped brackets, which enhance the safety of the structure. In the meantime, adjustable armrests, quick disassembly of footrest, positioning of Controller and other functions reflect the user-friendly design, and provide great convenience for the user.

Used in Patient Indoor-outdoor locomotion, Patient transfer.

Also known as Patient transfer wheelchair, Transport wheelchairs.

## HBD91-MW8 MANUAL WHEELCHAIR

It is a type of steady and simple Manual Wheelchair



### SPECIFICATIONS

Model	HBD91-MW8
Wheelchair length	1180 mm
Opened wheelchair width	660 mm
Wheelchair Height	1320 mm
Closed wheelchair width	360 mm
Seat Width	460 mm
Front wheel	8" inch
Rear wheel	24" inch
Seat Height	520 mm
Seat depth	460 mm
Backrest Height	490 mm
Armrest Height	230 mm
Distance between armrest	450 mm
Height between footrests and ground	115 mm
Net Weight	23.35 Kg





**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)