

# ELECTRIC WHEELCHAIR HBD92-EW2



# ELECTRIC WHEELCHAIR HBD92-EW2

Wheelchair is mainly used as a manually operated or power-driven device designed with double cross-shaped brackets, which enhance the safety of the structure. In the meantime, adjustable armrests, quick disassembly of footrest, positioning of Controller and other functions reflect the user-friendly design, and provide great convenience for the user.

Used in Seat elevation, Special operator controls for patient.

Also known as Motorized wheelchair, Electric-powered wheelchair.

## HBD92-EW2 ELECTRIC WHEELCHAIR

It is a type of steady, simple and Practical Electric Wheelchair



### SPECIFICATIONS

Model	HBD92-EW2
Total Length	1000 mm
Total Width	620 mm
Total Height	890 mm
Seat Depth	400 mm
Seat Width	460 mm
Seat Height	460 mm
Speed per Hour	> 6 Km/h
Maximum Load	160 Kg
Maximum Speed	7.2 Km/h
Mileage	17 Km
Damping System	2
Motor Power	24/250x2 V/W
Charger	24 (V2/A)
Km of Endurance	> 20 Km
Battery	24/10 ( V/Ah )
Charging Time	3-5 (h)
Front Wheel	8" inch
Rear Wheel	12" inch
Wheel Base	650 mm
Maximum Gradeability	12 °
Net Weight	25 Kg





**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)