

BULLET ICE MAKER ICM11-025



BULLET ICE MAKER ICM11-025

Compact machine with ice making capacity and innovative technology. It offers durability, efficiency and saving. Using quality assured components, it gives excellent functionality and low maintenance.

Used in Beverage industries, Biochemical and biological studies, Refrigerated perishables, Food industries, Medical laboratories, Healthcare, Pharmaceutical.

Also known as Ice Maker, Laboratory Bullet Ice Maker.

ICM11-025 BULLET ICE MAKER

Vertical stirring engine model

Refrigeration system produces uniform shape ice bullets

Bell-shaped ice with a hexagon hole

Highly durable stainless steel 304 2B built-up

Automatic controller increases working reliability

Large ice storage bin

Eco-friendly, energy saving non-Freon refrigerant



SPECIFICATIONS

Model	ICM11-025
Capacity	25 kg/24 h
Storage Bin Capacity	7 kg
Refrigerant	R134a
Condenser Type	Air Cooled
Overall Dimension	400x543x670 mm
Packing Dimension	470x615x730 mm
Weight	30/33 kg
Power	240 W
Power Supply	220V /110 V, 50 / 60 Hz

OPTIONAL ACCESSORIES

Accessory Code	Name	Description
LS54713	Water filter	



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com