

BULLET ICE MAKER ICM12-050



BULLET ICE MAKER ICM12-050

Compact machine with ice making capacity and innovative technology. It offers durability, efficiency and saving. Using quality assured components, it gives excellent functionality and low maintenance.

Used in Beverage industries, Biochemical and biological studies, Refrigerated perishables, Food industries, Medical laboratories, Healthcare, Pharmaceutical.

Also known as Ice Maker, Laboratory Bullet Ice Maker.

ICM12-050 BULLET ICE MAKER

Vertical stirring engine model

Refrigeration system produces uniform shape ice bullets

Bell-shaped ice with a hexagon hole

Highly durable stainless steel 304 2B built-up

Automatic controller increases working reliability

Large ice storage bin

Eco-friendly, energy saving non-Freon refrigerant



SPECIFICATIONS

Model	ICM12-050
Capacity	50 kg/24 h
Storage Bin Capacity	14 kg
Refrigerant	R134a
Condenser Type	Air Cooled
Overall Dimension	400x510x825 mm
Packing Dimension	470x590x875 mm
Weight	35/38 kg
Power	380 W
Power Supply	220V /110 V, 50 / 60 Hz

OPTIONAL ACCESSORIES

Accessory Code	Name	Description
LS54713	Water filter	



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com