

# BULLET ICE MAKER ICM12-120



# BULLET ICE MAKER ICM12-120

Compact machine with ice making capacity and innovative technology. It offers durability, efficiency and saving. Using quality assured components, it gives excellent functionality and low maintenance.

Used in Beverage industries, Biochemical and biological studies, Refrigerated perishables, Food industries, Medical laboratories, Healthcare, Pharmaceutical.

Also known as Ice Maker, Laboratory Bullet Ice Maker.

## ICM12-120 BULLET ICE MAKER

Transparent Bullet Ice

Air Cooling

Stainless Steel 304 Exterior

Adjustable Feet

Slide Away Access Door

World Famous Brand Compressor

Fully Automatic Controller

Improved Design For Higher Reliability



## SPECIFICATIONS

Model	ICM12-120
Capacity	120 kg/24 h
Storage Bin Capacity	55 kg
Refrigerant	R134a
Condenser Type	Air Cooled
Overall Dimension	680x611x1050 mm
Packing Dimension	750x675x1150 mm
Weight	70/75 kg
Power	800 W
Power Supply	220 V /110 V, 50 / 60 Hz



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)