

# CUBE ICE MAKER

## ICM21-0300



# CUBE ICE MAKER ICM21-0300

Ice-making speed is highly vigorous and it is energy efficient. Due to self contained storage bin, it beocmes more spacious and uniform cube shaped ice are produced rapidly.

Used in Healthcare, Food industries, Refrigerated perishables, Biochemical and biological studies, Medical laboratories, Beverage industries, Pharmaceutical.

Also known as Ice Cuber Machine, Laboratory Cube Ice Maker.

## ICM21-0300 CUBE ICE MAKER

High grade stainless steel SUS304 is used which is anticorrosive and durable. The compact structure is simple and space-saving.

Fully automatically computer controlled, convenient for use.

Automatic detection and protection in case of operation, water deficiency, ice full and malfunction.

Highly efficient fluoride-free compressor, energy saving and environment-friendly.

Suitable for Pubs, KTV and dancing halls, gyms, cold drinks shops, milk tea shops, cafes, convenience stores, restaurants and other entertainment venues.

Split-type structure, large ice bin, large ice storage capacity and strong applicability.



## SPECIFICATIONS

Model	ICM21-0300
Capacity	300 kg/24 h
Storage Bin Capacity	210 kg
Refrigerant	R134a / R404a
Condenser Type	Water cooled
Machine Head	760x600x560 mm
Cabinet	760x800x980 mm
Packing Machine Head	845x705x850 mm
Packing Cabinet	800x850x1060 mm
Weight	98/107 kg
Power Supply	220 V, 60 Hz



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)