

GEL CARD INCUBATOR INC31-24GC



GEL CARD INCUBATOR INC31-24GC

It can be suitable for many types of gel cards from different manufacturers. Feature like audible programmable alarm to notify end of incubation period makes it an ideal choice.

Used in Hospitals, Research, Laboratory.

Also known as Laboratory Gel Card Incubator.

INC31-24GC GEL CARD INCUBATOR

2 independent incubation zones with a capacity for 12 Gel cards

Incubator for 24-Gel Cards in racks at 37°C

Digital interface visor control and visualization of centrifuge speed and time remaining for processing



SPECIFICATIONS

Model	INC31-24GC
Temperature Range	RT +5-60°C
Temperature Accuracy	±1°C
Temperature Fluctuation	±0.5%
CO2 Range	0~20%
CO2 Recovery Time	≤concentration x 1.2 min
Jacket Type	Air jacket type
Overall Dimension (WxDxH)	220x380x140 mm
Power	300 W
Power Supply	AC220 V 3 A



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com