

GEL CARD INCUBATOR INC32-12GC



GEL CARD INCUBATOR INC32-12GC

It can be suitable for many types of gel cards from different manufacturers. Feature like audible programmable alarm to notify end of incubation period makes it an ideal choice.

Used in Hospitals, Research, Laboratory.

Also known as Laboratory Gel Card Incubator.

INC32-12GC GEL CARD INCUBATOR

Digital interface visor control and visualization of centrifuge speed and time remaining for processing

Automatic balance control

Pre-defined centrifugal parameters for optimum Gel cards processing; acceleration, speed, braking and time

Programmable audible alarms for end of centrifuge time periods



SPECIFICATIONS

Model	INC32-12GC
Capacity	12 Gel Cards
Speed Range	800~1800 rpm
RCF. Range	54~272 × g
Timer Range	1~20 min
Overall Dimension (WxDxH)	440 × 340 × 170 mm
Noise Level	≤55 db
Weight	7 kg
Power	85 W
Power Supply	AC100-240 V



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com