

GEL CARD INCUBATOR INC33-24GC



GEL CARD INCUBATOR INC33-24GC

It can be suitable for many types of gel cards from different manufacturers. Feature like audible programmable alarm to notify end of incubation period makes it an ideal choice.

Used in Hospitals, Research, Laboratory.

Also known as Laboratory Gel Card Incubator.

INC33-24GC GEL CARD INCUBATOR

Digital interface visor control and visualization of centrifuge speed and time remaining for processing

Automatic balance control

Pre-defined centrifugal parameters for optimum Gel cards processing; acceleration, speed, braking and time

Programmable audible alarms for end of centrifuge time periods



SPECIFICATIONS

Model	INC33-24GC
Capacity	24 Gel Cards
Speed Range	500~1500 rpm
RCF. Range	33~294 × g
Timer Range	1~20 min
Overall Dimension (WxDxH)	530 × 410 × 180 mm
Noise Level	≤55 db
Weight	11.5 kg
Power	85 W
Power Supply	AC100-240 V



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com