

BIOLOGICAL INDICATOR INCUBATOR INC52-080



BIOLOGICAL INDICATOR INCUBATOR INC52-080

Biological Indicator Incubator is supporting equipment for bio-indicator. It is small in size, less weight and is suitable for pressure steam sterilizer biological indicator, ethylene oxide sterilization biological indicator and other biological indicator sterilization, effect detection and cultivation.

Used in Environmental protection, Research institutes, University laboratories.

Also known as Laboratory Incubator.

INC52-080 BIOLOGICAL INDICATOR INCUBATOR

Designed for incubation of EO biological indicators

Customized wells for placement of biological indicator vials

Easy to set incubation temperature

Aluminium blocks available to meet different experimental requirements



SPECIFICATIONS

Model	INC52-080
Temperature Range	RT+5-80°C
Temperature Accuracy	0.35°C at 57°C
Temperature Uniformity	0.5°C
Temperature Resolution	0.1°C
Time for Heating	≤6 min (20°C to 57°C)
Timer Range	1 min-99h 59 min
Overall Dimension (WxDxH)	110x150x80 mm
Display	LCD
Weight	0.5/1 kg
Power	120 W
Power Supply	240 V, 60 Hz

OPTIONAL ACCESSORIES

Accessory Code	Name	Description
LS52989	Block C	6xΦ15 mm
LS53001	Block D	6xΦ16 mm
LS53013	Block E	24xΦ8.8 mm
LS53025	Block F	18xΦ9 mm
LS53049	Block B	15xΦ12 mm
LS54430	Block A	15 x Φ10 mm



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: info@biozef.com | Website: biozef.com