

# BIOLOGICAL INDICATOR INCUBATOR INC52-100



# BIOLOGICAL INDICATOR INCUBATOR INC52-100

Biological Indicator Incubator is supporting equipment for bio-indicator. It is small in size, less weight and is suitable for pressure steam sterilizer biological indicator, ethylene oxide sterilization biological indicator and other biological indicator sterilization, effect detection and cultivation.

Used in Environmental protection, Research institutes, University laboratories.

Also known as Laboratory Incubator.

## INC52-100 BIOLOGICAL INDICATOR INCUBATOR

Designed for incubation of EO biological indicators

Customized wells for placement of biological indicator vials

Easy to set incubation temperature

Aluminium blocks available to meet different experimental requirements



### SPECIFICATIONS

Model	INC52-100
Temperature Range	RT+5-100°C
Temperature Accuracy	0.35°C at 57°C
Temperature Uniformity	0.5°C
Temperature Resolution	0.1°C
Time for Heating	≤6 min (20°C to 57°C)
Timer Range	1 min-99h 59 min
Overall Dimension (WxDxH)	120x152x112 mm
Display	LCD
Weight	0.85/1.4 kg
Power	120 W
Power Supply	240 V, 60 Hz

### OPTIONAL ACCESSORIES

Accessory Code	Name	Description
LS52989	Block C	6xΦ15 mm
LS53001	Block D	6xΦ16 mm
LS53013	Block E	24xΦ8.8 mm
LS53025	Block F	18xΦ9 mm
LS53037	Block A	15xΦ10 mm
LS53049	Block B	15xΦ12 mm





**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)