

# BIOLOGICAL INDICATOR INCUBATOR INC53-057



# BIOLOGICAL INDICATOR INCUBATOR INC53-057

Biological Indicator Incubator is supporting equipment for bio-indicator. It is small in size, less weight and is suitable for pressure steam sterilizer biological indicator, ethylene oxide sterilization biological indicator and other biological indicator sterilization, effect detection and cultivation.

Used in Environmental protection, Research institutes, University laboratories.

Also known as Laboratory Incubator.

## INC53-057 BIOLOGICAL INDICATOR INCUBATOR

Compact, fits in the palm of your hand

LCD display. It is easy to setup and use

Different kinds of aluminum blocks are available

Clear cover ensures temperature uniformity

Portable, Optional 24V adapter



## SPECIFICATIONS

Model	INC53-057
Temperature Accuracy	0.35°C at 57°C
Temperature Control Accuracy	≤±0.5°C
Display Accuracy	0.1°C
Time for Heating	≤6 min (20°C to 57°C)
Timer Range	1 min-99 h 59 min
Overall Dimension (WxDxH)	120x152x112 mm

## OPTIONAL ACCESSORIES

Accessory Code	Name	Description
LS53257	Block	Ø10 mm X 15 tube Block
LS53269	Block	Ø12 mm X 12 tube Block
LS53281	Block	Ø15 mm X 6 tube Block
LS53293	Block	Ø16 mm X 16 tube Block
LS53305	Block	Ø9 mm X 24 tube Block





**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)