

# GENERAL PURPOSE INCUBATOR INC81-050



# GENERAL PURPOSE INCUBATOR INC81-050

General Purpose Incubator are designed to accommodate an ample range of laboratory applications and lends a contamination-free environment for cell and tissue cultures. These incubators serve the purpose for most biological analysis, research, general laboratory applications, biochemical studies, bacteriology, and hematological studies. Used in Quality Control, Medicine Processing Industries, Research, Pharmaceuticals, Microbial Culturing. Also known as Lab Incubator, Laboratory General Purpose Incubator.

## INC81-050 GENERAL PURPOSE INCUBATOR

- LCD PID control system
- International brand cooling fan
- Stainless steel chamber
- Cavity preheating technology and duct cycling system



## SPECIFICATIONS

Model	INC81-050
Capacity	50 L
Temperature Range	RT +5 ~ 65°C
Temperature Accuracy	±0.1°C
Temperature Fluctuation	±0.5%
Internal Dimension (WxDxH)	380x380x350 mm
Overall Dimension (WxDxH)	580x530x600 mm
Shelves	2
Weight	31 kg
Power Supply	220 V 50 Hz



### Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)