

GENERAL PURPOSE INCUBATOR INC81-080



GENERAL PURPOSE INCUBATOR INC81-080

General Purpose Incubator are designed to accommodate an ample range of laboratory applications and lends a contamination-free environment for cell and tissue cultures. These incubators serve the purpose for most biological analysis, research, general laboratory applications, biochemical studies, bacteriology, and hematological studies. Used in Quality Control, Medicine Processing Industries, Research, Pharmaceuticals, Microbial Culturing. Also known as Lab Incubator, Laboratory General Purpose Incubator.

INC81-080 GENERAL PURPOSE INCUBATOR

LCD PID control system
International brand cooling fan
Stainless steel chamber
Cavity preheating technology and duct cycling system



SPECIFICATIONS

| | |
|----------------------------|----------------|
| Model | INC81-080 |
| Capacity | 80 L |
| Temperature Range | RT +5 ~ 65°C |
| Temperature Accuracy | ±0.1°C |
| Temperature Fluctuation | ±0.5% |
| Internal Dimention (WxDxH) | 480x480x350 mm |
| Overall Dimension (WxDxH) | 680x630x600 mm |
| Shelves | 2 |
| Weight | 42 kg |
| Power Supply | 220 V 50 Hz |



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com