

# GENERAL PURPOSE INCUBATOR INC81-160



# GENERAL PURPOSE INCUBATOR INC81-160

General Purpose Incubator are designed to accommodate an ample range of laboratory applications and lends a contamination-free environment for cell and tissue cultures. These incubators serve the purpose for most biological analysis, research, general laboratory applications, biochemical studies, bacteriology, and hematological studies. Used in Quality Control, Medicine Processing Industries, Research, Pharmaceuticals, Microbial Culturing. Also known as Lab Incubator, Laboratory General Purpose Incubator.

## INC81-160 GENERAL PURPOSE INCUBATOR

- LCD PID control system
- International brand cooling fan
- Stainless steel chamber
- Cavity preheating technology and duct cycling system



### SPECIFICATIONS

Model	INC81-160
Capacity	160 L
Temperature Range	RT +5 ~ 65°C
Temperature Accuracy	±0.1°C
Temperature Fluctuation	±0.5%
Internal Dimension (WxDxH)	600x600x450 mm
Overall Dimension (WxDxH)	800x750x700 mm
Shelves	2
Weight	57 kg
Power Supply	220 V 50 Hz



### Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)