

# GENERAL PURPOSE INCUBATOR INC82-100



# GENERAL PURPOSE INCUBATOR INC82-100

General Purpose Incubator are designed to accommodate an ample range of laboratory applications and lends a contamination-free environment for cell and tissue cultures. These incubators serve the purpose for most biological analysis, research, general laboratory applications, biochemical studies, bacteriology, and hematological studies. Used in Quality Control, Medicine Processing Industries, Research, Pharmaceuticals, Microbial Culturing. Also known as Lab Incubator, Laboratory General Purpose Incubator.

## INC82-100 GENERAL PURPOSE INCUBATOR

- Balance Refrigeration
- Adopt without fluorine R134a refrigerant technology
- Duct circulation
- Stainless steel inner
- Programmable intelligent PID control
- Usd large capacity humidifier
- Double door design
- Lighting device
- With UV lamp



### SPECIFICATIONS

Model	INC82-100
Capacity	100 L
Temperature Range	0-60°C
Temperature Accuracy	±0.1°C
Temperature Fluctuation	±1%
Humidity Range	50%-90 %
Humidity Fluctuation	±5%
Timer Range	0-999 min
Internal Dimension (WxDxH)	450×380x580 mm
Overall Dimension (WxDxH)	590×620x1100 mm
Package Size	670x700x1220 mm
Weight	93/100 kg
Power	400 W
Power Supply	220 V



### Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)