

GENERAL PURPOSE INCUBATOR INC82-150



GENERAL PURPOSE INCUBATOR INC82-150

General Purpose Incubator are designed to accommodate an ample range of laboratory applications and lends a contamination-free environment for cell and tissue cultures. These incubators serve the purpose for most biological analysis, research, general laboratory applications, biochemical studies, bacteriology, and hematological studies. Used in Quality Control, Medicine Processing Industries, Research, Pharmaceuticals, Microbial Culturing. Also known as Lab Incubator, Laboratory General Purpose Incubator.

INC82-150 GENERAL PURPOSE INCUBATOR

Balance Refrigeration

Adopt without fluorine R134a refrigerant technology

Duct circulation

Stainless steel inner

Programmable intelligent PID control

Usd large capacity humidifier

Double door design

Lighting device

With UV lamp



SPECIFICATIONS

Model	INC82-150
Capacity	150 L
Temperature Range	0-60°C
Temperature Accuracy	±0.1°C
Temperature Fluctuation	±1%
Humidity Range	50%-90 %
Humidity Fluctuation	±5%
Timer Range	0-999 min
Internal Dimention (WxDxH)	500×380x800 mm
Overall Dimension (WxDxH)	640×660x1350 mm
Package Size	720x740x1470 mm
Weight	105/115 kg
Power	600 W
Power Supply	220 V



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com