

GENERAL PURPOSE INCUBATOR INC82-250



GENERAL PURPOSE INCUBATOR INC82-250

General Purpose Incubator are designed to accommodate an ample range of laboratory applications and lends a contamination-free environment for cell and tissue cultures. These incubators serve the purpose for most biological analysis, research, general laboratory applications, biochemical studies, bacteriology, and hematological studies. Used in Quality Control, Medicine Processing Industries, Research, Pharmaceuticals, Microbial Culturing. Also known as Lab Incubator, Laboratory General Purpose Incubator.

INC82-250 GENERAL PURPOSE INCUBATOR

- Balance Refrigeration
- Adopt without fluorine R134a refrigrant technology
- Duct circulation
- Stainless steel inner
- Programmable intelligent PID control
- Usd large capacity humidifier
- Double door design
- Lighting device
- With UV lamp



SPECIFICATIONS

Model	INC82-250
Capacity	250 L
Temperature Range	0-60°C
Temperature Accuracy	±0.1°C
Temperature Fluctuation	±1%
Humidity Range	50%-90 %
Humidity Fluctuation	±5%
Timer Range	0-999 min
Internal Dimention (WxDxH)	550×450x1000 mm
Overall Dimension (WxDxH)	690×660x1550 mm
Package Size	770x740x1670 mm
Weight	133/145 kg
Power	800 W
Power Supply	220 V



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com