

# AIR-JACKETED MULTI GAS INCUBATOR INCB1-151L



# AIR-JACKETED MULTI GAS INCUBATOR

## INCB1-151L

Designed for a broad range of applications in the biomedical, pharmaceutical, medical and clinical research. Air-jacketed Multi gas Incubator system help to facilitate a stable cell culture environment. Easy to clean and minimizes recovery time and the risk of contamination.

Used in Biomedical, Pharmaceutical, Medical, Clinical research, Laboratory.

Also known as Laboratory Air-jacketed Multi gas Incubator, Multi gas Incubator Air-jacketed.

## INCB1-151L AIR-JACKETED MULTI GAS INCUBATOR

151 liters internal volume

Zirconium oxide Sensor

IR CO<sub>2</sub> sensor

Water Level Alarm

Integrated electric siphon pump

Auto-zero

Temperature deviation (at 37° C)

Moist-heat disinfection at 90°C

Relative humidity ≥95%



## SPECIFICATIONS

Model	INCB1-151L
Capacity	151 L
Temperature Sensors	PT1000
CO <sub>2</sub> Control Method	Microprocessor PID
Shelve Dimension	423×445 mm
Internal Dimension (WxDxH)	470×530×607 mm
Overall Dimension (WxDxH)	637×768×869 mm
Shelves	3 standard, 10 maximum
Weight	80 kg
Max. Load per Shelf	10 kg
Power	≤650 VA+10%
Power Supply	220 V±10%/ 50 Hz

## ACCESSORIES

Accessory Code	Name	Description	3 Sets	1 Set
LS53974	Standard shelf		3 Sets	
LS53998	CO <sub>2</sub> steel cylinder auto-switch device			1 Set

## OPTIONAL ACCESSORIES

Accessory Code	Name	Description
LS53650	Air connector	
LS53662	Shelf support	
LS53686	Small HEPA filter	
LS53758	CO2 regulator	

Accessory Code	Name	Description
LS53782	O2 regulator	
LS53794	Shelf	
LS53806	N2 regulator	



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)