

PLATELET INCUBATOR INCC1-060



PLATELET INCUBATOR INCC1-060

Platelet Incubator is designed to store platelet concentrates in continuous motion at a controlled temperature. It provides precise and steady storage conditions for maintaining the viability of platelet concentrates. Used in Biochemical Studies, Hematological Studies, Bacterial Culturing, Research, Microbiological Determinations, Pharmaceutical Stability Assays, Food Processing Quality Control. Also known as Laboratory Platelet Incubator.

INCC1-060 PLATELET INCUBATOR

Microcomputer temperature controller

LCD Display

'Set temperature' and 'Actual temperature' display on double screen, clear and straight forward.

Air forced convection design, ensures the temperature uniform in the chamber.

Over temperature / Low temperature alarm, abnormal shutdown power outage alarm.

Sets in exchange stomatal, to ensure that the air is fresh in the cabinet.

UV Lamp in the cabinet.

Automatic defrost function.

Big observation window, working condition will be clear at a glance.

Stainless steel inner chamber material, easy to clean.

Heat preservation material, uses polyurethane foam.

Closed type compressor.

Less noise.



SPECIFICATIONS

Model	INCC1-060
Capacity	10/15/20 floors (40-60bags)
Temperature Range	15-35°C
Temperature Accuracy	±0.1
Timer Range	0-999
Variation (%)	±0.5
Internal Dimension (WxDxH)	590x390x820 mm
Overall Dimension (WxDxH)	670x600x1500 mm
Tray Size (WxDxH)	480x300x680 mm
Weight	180 Kg
Power	0.8 KW



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com