

PH METER MET21-14



PH METER MET21-14

The pH scale typically ranges from 0 to 14. According to the Nernst equation pH values are affected by temperature, hence the temperature compensation is a vital function of any pH Meter to ensure accurate pH measurement.

MET21-14 PH METER



Product Image Coming Soon

SPECIFICATIONS

Model	MET21-14
pH	
Range	-2.00 to 14.00 pH
Resolution	0.01 pH
Accuracy	± 0.01 pH
Calibration Points	Up to 2
Standard Recognition	NIST buffers
mV	
Range	-1800 to 1800 mV
Resolution	1 mV
Accuracy	± 0.1 %FS
Conductivity	
Range	/
Resolution	/
Accuracy	/
Reference Temperature	/
Measurement	
Reading Mode	Continuous
Reading Prompts	Reading, Stable
Temp. Compensation	MTC
Inputs	
pH Electrode	BNC(Q9)
Conductivity Probe	/
Display Options	
Backlight	Yes
Auto Shutdown	300, 600, 1200, 1800, 3600 sec., off

IP Rating	IP54
General	
Power	AC Adapter, 100-240 V AC input, DC 9V output
Dimensions	200×160×63 mm
Weight	600 g (1.32 lb)



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: info@biozef.com | Website: biozef.com