

# TUBE ROLLER MIX15-80T



# TUBE ROLLER MIX15-80T

DC brushless motor provides reliable and quiet running, maintenance free and can be used in cold rooms or in incubators, Majorly used in Industrial, Clinical Laboratories, Academic Facilities, Environmental Applications, Food Technology and Wastewater Facilities.

Used in Industrial, Clinical Laboratories, Academic Facilities, Environmental Applications, Food Technology, Wastewater Facilities.

Also known as Tube Rotator, Rotator, Laboratory Rotator.

## MIX15-80T TUBE ROLLER

End-over-end or rolling action, gives gentle but effective mixing

Adjustable speed range of 0 to 80 rpm

Used with a wide variety of accessories

Continuous quiet operation



## SPECIFICATIONS

Model	MIX15-80T
Speed Range	0-80 rpm
Permissible Ambient Temperature	5-40°C
Permissible Relative Humidity	0.8
Tilt Angle (°)	90
Operation	Continuous
Overall Dimension	300Dx220Wx310H mm
Motor Type	DC motor
Protection Class	IP21
Weight	2 kg
Power	20 W
Power Supply	100-240V, 50/60Hz

## OPTIONAL ACCESSORIES

Accessory Code	Name	Description
LS52589	Disk and Rotisserie accessory	For 1.5ml x 12 centrifuge tubes
LS52602	Disk and Rotisserie accessory	For 15ml x 12 centrifuge tubes
LS52615	Disk and Rotisserie accessory	For 50ml x 12 centrifuge tubes





**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)