

# VORTEX MIXER MIX31-28



# VORTEX MIXER MIX31-28

Constructed of rugged, corrosion-resistant materials, the instruments endure heavy everyday usage. The ergonomic low profile design makes everyday vortexing comfortable for the user. Unnecessary movement has been avoided during use due to robust die-cast body.

Used in Mixing, Laboratory, Research, Medical.

Also known as Laboratory Vortex Mixer.

## MIX31-28 VORTEX MIXER

The maximum diameter of the mixed tube is up to 30mm

Suit for mixing tubes / centrifugal tubes with prominent mixing effect

Outer shell and touch concave are made of plastic which is resistance to chemical attack

The thick base makes operation very steady and reliable

12V low voltage power adapter built in, very smart and portable



## SPECIFICATIONS

Model	MIX31-28
Speed Range	2800 rpm
Capacity	up to 50 ml (1 tube)
Orbital Diameter	4.5 mm
Motion	Circle Shaking (Vortexing)
1.2	Motor rating output [W&#93
Overall Dimension	110 x 70 mm
Weight	1.2 kg
Power	15 W
Power Supply	AC 90~230V



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)