

VORTEX MIXER MIX31-28A



VORTEX MIXER MIX31-28A

Constructed of rugged, corrosion-resistant materials, the instruments endure heavy everyday usage. The ergonomic low profile design makes everyday vortexing comfortable for the user. Unnecessary movement has been avoided during use due to robust die-cast body.

Used in Mixing, Laboratory, Research, Medical.

Also known as Laboratory Vortex Mixer.

MIX31-28A VORTEX MIXER

Touch and continuous operations, stepless speed regulation

Used for various mixing applications with optional adapters, handles up to 48 tubes

Specially designed vacuum suction feet hold to the work surface and prevent the unit from movement

Aluminium cast base avoids unnecessary "walking"



SPECIFICATIONS

Model	MIX31-28A
Speed Range	0-2800 rpm
Capacity	up to 50 ml (1 tube)
Orbital Diameter	4.5 mm
Motion	Circle Shaking (Vortexing)
1.2	Motor rating output [W]
Overall Dimension	110 x 70 mm
Weight	1.2 kg
Power	15 W
Power Supply	AC 90~230V



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com