

VORTEX MIXER MIX37-25A



VORTEX MIXER MIX37-25A

Constructed of rugged, corrosion-resistant materials, the instruments endure heavy everyday usage. The ergonomic low profile design makes everyday vortexing comfortable for the user. Unnecessary movement has been avoided during use due to robust die-cast body.

Used in Mixing, Laboratory, Research, Medical.

Also known as Laboratory Vortex Mixer.

MIX37-25A VORTEX MIXER

Suit for mixing tubes / centrifugal tubes with prominent mixing effect

Outer shell and touch concave are made of plastic which is resistance to chemical attack

The thick base makes operation very steady and reliable

12V low voltage power adapter built in, very smart and portable

Can adjust speed to meet different experimental requests



SPECIFICATIONS

| | |
|---------------------------------|-----------------------------|
| Model | MIX37-25A |
| Speed Range | 0-2500 rpm |
| Permissible Ambient Temperature | 5-40°C |
| Permissible Relative Humidity | 0.8 |
| Orbital Diameter | 4 mm |
| Motion | Orbital |
| 58 | Motor rating output [W] |
| Operation | Continuous/ touch operation |
| Overall Dimension | 127Wx130Hx160D mm |
| Speed Display | Scale |
| Protection Class | IP21 |
| Weight | 3.5 kg |
| Power | 60 W |
| Power Supply | 110-120V /220-230V, 50/60Hz |

OPTIONAL ACCESSORIES

| Accessory Code | Name | Description |
|----------------|--|--|
| LS52332 | Rubber vacuum suction foot with screw, 1pc | |
| LS52371 | Platform pad | For $\varnothing 99\text{mm}$ tubes and small vessels |

| Accessory Code | Name | Description |
|----------------|--------------|---|
| LS52462 | Tube adapter | For 12 holes test tubes, $\varnothing 12\text{mm}$, used with MX1027/ MX1026/PS1.1 |

| | | |
|---------|------------------|--|
| LS52384 | Standard top cup | For \varnothing30mm tubes and small vessels |
| LS52436 | Tube adapter | For 48 holes test tubes, \varnothing 6mm, used with MX1027/MX1026/PS1.1 |
| LS52449 | Tube adapter | For 15 holes test tubes, \varnothing 10mm, used with MX1027/MX1026/PS1.1 |

| | | |
|---------|---------------------|--|
| LS52475 | Tube adapter | For 8 holes test tubes, \varnothing 16mm, used with MX1027/MX1026/PS1.1 |
| LS52488 | Tube adapter | For 8 holes test tubes, \varnothing 20mm, used with MX1027/MX1026/PS1.1 |
| LS52501 | Universal top plate | 100mm plate attachment , free of Rubber vacuum suction foot (18900158, 4pcs) |
| LS52527 | Tube holding rod | Used with MX1028 to MX1032 |



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com