

# BIOLOGICAL MICROSCOPE MSC21-13



# BIOLOGICAL MICROSCOPE MSC21-13

Digital Biological Microscope is an optical microscope that is primarily designed to observe cells, tissues as well as other biological specimens. Magnification range anywhere from 10x – 1,000x or more by attaching multiple objective lenses.

Used in University, Medical, Laboratory, Research, Pharmaceutical.

## MSC21-13 BIOLOGICAL MICROSCOPE

The integration of structure design to save space.

Valid mould proof design.

It can be used under the high temperature and high humidity environment.



## SPECIFICATIONS

Model	MSC21-13
Viewing Head	Binocular head
Fluorescent Straight Illumination	B(Blue), G(Green) exciting light filter system O ordinary light system 100 W HBO super high-pressured spherical mercury lamp
Protection Screen	Screen to resist the UV
Power	AC input 220 V
Connecting Wire	With 2 connectors
Fluorescent Free Objective	10X, 20X, 40X (S), 100X (S, Oil)
Immersion Oil	Fluorescent free oil
Power Supply	AC110 / 220 V ± 10 % , 50/60 Hz
Package Size	440x340x220 mm
Gross Weight	17 kg



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)