

BIOLOGICAL MICROSCOPE MSC21-21



BIOLOGICAL MICROSCOPE MSC21-21

Digital Biological Microscope is an optical microscope that is primarily designed to observe cells, tissues as well as other biological specimens. Magnification range anywhere from 10x – 1,000x or more by attaching multiple objective lenses.

Used in University, Medical, Laboratory, Research, Pharmaceutical.

MSC21-21 BIOLOGICAL MICROSCOPE

The integration of structure design to save space.

Valid mould proof design.

It can be used under the high temperature and high humidity environment.



Product Image Coming Soon

SPECIFICATIONS

Model	MSC21-21
Viewing Head	Monocular Head Inclined at 45°, 360° rotatable
Achromatic Objectives	4X, 10X, 40X (Spring) , 100X (Spring, Oil)
Eyepiece	Wide Field Eyepiece:WF10 X, (WF16 X optional)
Stage	Double Layer Mechanical Stage Size 120x125 mm
Focusing	Coaxial Coarse and Fine Adjustment, Focusing Range 30 mm, Focusing Interval 0.002 mm
Condenser	Abbe NA=1.25 with Iris Diaphragm&Filter
Illumination	LED(1 W) Lamp 220 V or 110 V, (LED 3 W lamp optional)
Power Supply	AC110/220 V ± 10 %, 50/60 Hz
Package Size	450x360x240 mm
Gross Weight	5 kg



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com