

BIOLOGICAL MICROSCOPE MSC21-22



BIOLOGICAL MICROSCOPE MSC21-22

Digital Biological Microscope is an optical microscope that is primarily designed to observe cells, tissues as well as other biological specimens. Magnification range anywhere from 10x – 1,000x or more by attaching multiple objective lenses.

Used in University, Medical, Laboratory, Research, Pharmaceutical.

MSC21-22 BIOLOGICAL MICROSCOPE

The integration of structure design to save space.

Valid mould proof design.

It can be used under the high temperature and high humidity environment.



SPECIFICATIONS

Model	MSC21-22
Viewing Head	Compensation Free Binocular Head Inclined at 30°, 360° rotatable
Achromatic Objectives	4X, 10X, 40X (Spring) , 100X (Spring, Oil)
Eyepiece	Wide Field Eyepiece:WF10 X, (WF16 X optional)
Stage	Double Layer Mechanical Stage Size 120x125 mm
Focusing	Coaxial Coarse and Fine Adjustment, Focusing Range 30 mm, Focusing Interval 0.002 mm
Condenser	Abbe NA=1.25 with Iris Diaphragm&Filter
Illumination	LED(1 W) Lamp 220 V or 110 V, (LED 3 W lamp optional)
Power Supply	AC110/220 V \pm 10 %, 50/60 Hz
Package Size	450x360x240 mm
Gross Weight	5 kg



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com