

STEREO ZOOM MICROSCOPE MSC41-040ST



STEREO ZOOM MICROSCOPE MSC41-040ST

Stereo Zoom Microscope also known as the dissecting microscope is type of optical microscope that allows the user to see a three-dimensional view of a specimen.

Used in Pharmaceutical, University, Medical, Laboratory, Research.

MSC41-040ST STEREO ZOOM MICROSCOPE

Excellent optical function with infinite optical system.

Innovative stand structure and sharp image display.

Convenient operation and special design for viewing incubating cell tissue.

Valid mould proof design.

It can be used under the high temperature and high humidity environment.



SPECIFICATIONS

Model	MSC41-040ST
Illumination	Tranilluminator / inclined light
View Head And Body	Inclined viewing Up-rotatable body
Objective	2X&4X, 90°rotated
Working Distance	57 mm
Eyepiece	WF10X/20
Power Supply	AC110/220 V ± 10%, 50/60 Hz
Package Size	360x170x290 mm
Gross Weight	6 kg



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com