

STEREO ZOOM MICROSCOPE MSC41-045MT



STEREO ZOOM MICROSCOPE MSC41-045MT

Stereo Zoom Microscope also known as the dissecting microscope is type of optical microscope that allows the user to see a three-dimensional view of a specimen.

Used in Pharmaceutical, University, Medical, Laboratory, Research.

MSC41-045MT STEREO ZOOM MICROSCOPE

Excellent optical function with infinite optical system.

Innovative stand structure and sharp image display.

Convenient operation and special design for viewing incubating cell tissue.

Valid mould proof design.

It can be used under the high temperature and high humidity environment.



SPECIFICATIONS

| Model | MSC41-045MT |
|--|-------------|
| Viewing Head | |
| Binocular view head inclined at 60° | No |
| Trinocular view head inclined at 45° | Yes |
| Digital view head inclined at 45° | Yes |
| Eyepiece | |
| High point wide field adjustable eyepiece WF10x22 mm | Yes |
| Zoom Objective | |
| 0.67 X—4.5 X | Yes |
| Zoom Ratio | |
| 6.7:1 | Yes |
| Auxiliary Objective | |
| 0.5X /177 mm, 1.5X /47 mm, 2X/26 mm | Yes |
| Working Distance | |
| 115 mm | Yes |
| Focusing Bracket | |
| Focusing frame and adjustment knob | Yes |
| Interpupillary Distance | |
| 55-75 mm / 52-75 mm | Yes |
| Up Illumination | |
| LED lamp | Yes |
| Down Illumination | |
| LED lamp | Yes |
| Power Supply | |
| AC110/220 V ± 10%, 50/60 Hz | Yes |

| | |
|-------------------|----------------|
| Package | |
| Cartoon with foam | Yes |
| Package Size | 500x260x370 mm |
| Gross Weight | 8 kg |



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: info@biozef.com | Website: biozef.com