

# STEREO ZOOM MICROSCOPE MSC41-060ST



# STEREO ZOOM MICROSCOPE MSC41-060ST

Stereo Zoom Microscope also known as the dissecting microscope is type of optical microscope that allows the user to see a three-dimensional view of a specimen.

Used in Pharmaceutical, University, Medical, Laboratory, Research.

## MSC41-060ST STEREO ZOOM MICROSCOPE

Excellent optical function with infinite optical system.

Innovative stand structure and sharp image display.

Convenient operation and special design for viewing incubating cell tissue.

Valid mould proof design.

It can be used under the high temperature and high humidity environment.



## SPECIFICATIONS

Model	MSC41-060ST
Illumination	Natural light
View Head And Body	Inclined viewing Body 360°rotatable
Objective	2X&4X, 360°rotated
Working Distance	100 mm
Eyepiece	WF10X/20
Power Supply	AC110/220 V ± 10%, 50/60 Hz
Package Size	400x250x360 mm
Gross Weight	8 kg



**Biozef**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [info@biozef.com](mailto:info@biozef.com) | Website: [biozef.com](http://biozef.com)