

SINGLE ARM LED OPERATING LAMP OPL16-14



SINGLE ARM LED OPERATING LAMP OPL16-14

Challenging tasks require special solutions. The innovative Labstac LED operation light give you benefits of cutting-edge LED technology in a uniquely effective manner. We focus on providing LED lights that are more color-accurate, more flexible, more economical. We offer innovative light and system solutions that can be retrofitted easily and without special preparation.

Used in Dental implant surgery, Medical unit, Surgery, Operation theatre.

Also known as Surgical LED Light, LED operation theatre light, LED OT light.

OPL16-14 SINGLE ARM LED OPERATING LAMP

Streamlined, ultra-thin, enclosed head with reduced air resistance to minimise contamination

Digital controlled light intensity without grades of brightness and automatic memory for added flexibility

Optionally built in power supply systems avoiding complicated installation of wall mounted control boxes

Auxiliary light source included for laparoscopic operations

Rigid suspension arm with easy mobility and stable positioning

No infrared or ultraviolet radiation

Backlight compensation function

Endoscopy light source

Stepless adjusting light source



SPECIFICATIONS

Model	OPL16-14
Illumination	60000~120000 Lux
Diameter of spot	10-28 cm
Depth of light beam	≥700 mm
Color rendering Index (Ra)	100≥Ra≥85
Color Temperature (K)	4500±500 K



Biozef

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: info@biozef.com | Website: biozef.com